Bravo™

BLOOD GLUCOSE TEST STRIPS

IMPORTANT

Please read this information and the **Bravo™** Owner's booklet before using the **Bravo™** Test Strips.

INTENDED USE AND TEST PRINCIPLE

The **Bravo™** Test Strip works with the **Bravo™** Blood Glucose Meter to quantitatively measure glucose in whole blood. The **Bravo™** Blood Glucose Monitoring System is for self-testing outside the body (*in vitro* diagnostic use). The **Bravo™** Blood Glucose Monitoring System should not be used for the diagnosis of diabetes or for testing newborns. The glucose in a blood sample mixes with special test strip chemicals where a small electrical current is produced. The **Bravo™** Meter detects this electrical current and measures the amount of glucose in the blood sample.

STORAGE AND HANDLING

- Store test strips only in their original vial.
- Store vial in a cool and dry place between 1-30° C (34-86° F).
- Keep the vial away from direct sunlight or heat.
- Do not freeze test strips or vial.
- Close vial cap tightly immediately after removing a test strip.
- Use test strip immediately after taking it out of the vial.
- Handle test strips only with clean and dry hands.
- Do not bend, cut, or alter the test strips in any way.
- Do not force a test strip into the meter. Gently push it into the meter's test strip port.
- Apply *only* fresh whole capillary blood to the test strip.
- Do not use test strips past expiration or discard date. Use test strips within six (6) months of opening the vial. Record the discard date (the date the vial was opened plus six (6) months) on the vial label.

WARNINGS AND PRECAUTIONS

- Keep test strips and the test strip vial away from children. The test strips and vial cap may be choking hazards. Drying agents in the vial cap may be harmful if inhaled or swallowed or may cause skin or eye irritation.
- Test strips are for single use only. Do not reuse.
- If the test strip does not absorb the blood sample properly, please contact your authorized **Endomedical™** sales representative.

BLOOD SAMPLE COLLECTION PROCEDURE

Wash hands and sample site with soap and warm water. Rinse and dry thoroughly before collecting the blood sample with a lancing device.

Fingertip Site Blood Sampling

Unscrew lancing device tip. Place cocked lancing device against fingertip side. Press release button. Massage fingertip to obtain round drop (at least 0.5 μ L) of blood. Apply test strip tip to blood sample.

Alternative Site Blood Sampling (forearm, palm, thigh, calf)

Select a clean, soft and fleshy sample site area free of visible veins and hair and away from bones. Gently massage the sample site to help blood circulation to minimize result differences between fingertip and alternative site sampling. Firmly press and hold the lancing device against site. Wait until the skin surface under the lancing device changes color. Then press the release button while continuing to apply pressure. Keep holding the lancing device against your skin until sufficient (at least 0.5 μ L) blood is drawn. Carefully lift the lancing device away from your skin.

Acceptable situations for AST

- Fasting period
- Before a meal


Situations requiring fingertip test


- When the glucose levels are rapidly changing (during the two hours after a meal, insulin dose or exercise)
- When sick or when glucose levels seem quite lower than test value
- When hypoglycemia is not well recognized (hypoglycemic unawareness)
- When insulin has the biggest effect
- 2 hours after an insulin injection
- If you think your blood glucose is low
- If you think you will be driving a car or operating machinery

CAUTION

Alternative site and fingertip results may differ significantly due to rapid changes in the glucose level after meals or exercise, hypoglycemic symptoms, or effects of drugs such as insulin. Rubbing the test site before collecting blood reduces these differences. Use a fingertip sample site if you suffer from hypoglycemia or have experienced hypoglycemic shock or symptoms.

TEST PROCEDURE

- 1) Wash hands and sample site with soap and warm water. Rinse and dry thoroughly.
- 2) Insert the test strip into the port with contact bars facing upwards. Push the strip in gently until the meter beeps.
- 3) The  symbol will appear.

- 4) Use lancing device to get blood sample. Sample must be at least 0.5 μ L to fill the test strip confirmation window. When the  symbol appears on display, apply blood sample to edge of the narrow end of the test strip until meter beeps. If confirmation window has not filled, Er4 message or inaccurate test results will appear.
- 5) Meter will count down from five-to one (5 -to-1) on the display. Test result, time, and date will appear and automatically store in meter memory. Remove used test strip from port. Meter will turn off after three (3) seconds.

TEST RESULTS

The **Bravo™** Blood Glucose Meter will display results between 1.1-33.3 mmol/L.

Blood Glucose Targets for Most People with Diabetes

A recommended target for fasting blood glucose or glucose before meals is 4,0-7,0 mmol/L. A recommended target for blood glucose 2 hours after eating is 5,0 to 10,0 mmol/L.¹

Low Blood Glucose Results

If the test result is below 1.1 mmol/L, **Lo** will appear on the display indicating hypoglycemia (low blood glucose). You should follow the appropriate treatment recommendations of your healthcare professional.

High Blood Glucose Results

If the test result is above 33.3 mmol/L, **Hi** will appear on the display to indicate hyperglycemia (high blood glucose). If so, follow hyperglycemia treatment recommendations of your healthcare professional.

Unexpected Results

Low or high blood glucose readings can indicate a potentially serious medical condition. If your results are unusually high or low, or do not match the way you feel, repeat the test with a new test strip. If your reading is inconsistent with your symptoms or your result is less than 3.3 mmol/L or higher than 13.3 mmol/L, contact your healthcare professional.

Please note that:

- An abnormally high or low red blood cell count (hematocrit level over 60% or below 20%) may produce inaccurate results.
- Severe dehydration (excessive water loss) may cause false low results. If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- There is no effect from altitude up to 3,000 m (10,000 ft) above sea level.
- Interferences: Acetaminophen, ascorbic acid (vitamin C), uric acid and other reducing substances (when occurring in normal blood or normal therapeutic concentrations) do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurate high results.
- Blood samples that contain a high concentration of dissolved oxygen may lower the test result.
- Discard used test strips properly in a container for sharp objects.

METER AND TEST STRIP PERFORMANCE CHECK

The **Bravo™** Control Solution (Control A and/or B) contains a known amount of glucose that reacts with the **Bravo™** Test Strip in combination with the **Bravo™** Meters to make sure they are working properly together and the correct testing procedure is being followed.

You may run a check when you:

- Want to practice the test procedure using the control solution instead of blood.
- Use the meter for the first time.
- Open a new vial of test strips.
- Have symptoms that are inconsistent with your blood glucose test results.
- Believe your test results are not accurate.
- Suspect your meter and test strips are not performing properly.

If your control solution test results do not fall within the range printed on the test strip vial, repeat the test. Out of range results may be due to one or more of the following factors:

- Error in performing the test.
- Expired or contaminated control solution.
- Expired or damaged test strip.
- Failure to shake control solution bottle.
- Failure to discard first drop of control solution and wipe bottle tip clean.

If results continue to fall outside the range printed on the test strip vial, the **Bravo™** Test Strip and Meter may not be working properly. If so, do not use your system and contact your authorized **Endomedical™** sales representative.

CHEMICAL COMPOSITION

Each **Bravo™** Test Strips contains the following reagents:
Enzyme: Glucose oxidase (*Aspergillus* sp.): 2.7 units
Mediator: Hexaamineruthenium (III) chloride : 45.7 ug
Other ingredients: 1.6 ug

PERFORMANCE CHARACTERISTICS

The performance of **Bravo™** BGM System has been evaluated in laboratory and in clinical tests.

ACCURACY

Bravo™ BGM System is calibrated to yield results equivalent to plasma glucose concentrations. The accuracy of the **Bravo™** BGM System (Model GM505MA, GM505MB, GM505MC, GM505MD, GM505ME) was assessed by comparing blood glucose results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

Slope	0.97
Y-intercept	-0.133 mmol/L
Correlation coefficient (r)	0.993
Number of Sample	114
Range tested	2.5-32.5 mmol/L

Accuracy results for glucose concentration < 4.2 mmol/L		
Within± 0.28 mmol/L	Within± 0.56 mmol/L	Within ± 0.83 mmol/L
8/8 (100%)	8/8 (100%)	8/8 (100%)

Accuracy results for glucose concentration ≥ 4.2 mmol/L			
Within± 5%	Within± 10%	Within± 15%	Within± 20%
49/106 (46%)	82/106 (77%)	101/106 (95%)	106/106 (100%)

PRECISION

Precision studies were performed in a laboratory using the **Bravo™** BGM System.

Within Run Precision

*Blood avg.	2.8 mmol/L	SD=0.1 mmol/L
*Blood avg.	5.1 mmol/L	SD=0.2 mmol/L
*Blood avg.	7.8 mmol/L	CV=2.4%
*Blood avg.	11.3 mmol/L	CV=2.9%
*Blood avg.	19.2 mmol/L	CV=2.9%

Total Precision

*Control avg.	2.6 mmol/L	SD=0.1 mmol/L
*Control avg.	7.6 mmol/L	CV=3.1%
*Control avg.	18.7 mmol/L	CV=3.5%


This study shows that there could be variation of up to 3.5%.

Reference

1. Canadian Diabetes Association 2008 Clinical Practice Guidelines for the Prevention and Management of Diabetes in Canada.

These targets are intended to be a guide. Please discuss your blood glucose targets with your physician.

DESCRIPTION OF SYMBOL

	Consult instruction for use
	Temporary limitations
	<i>In vitro</i> diagnostic medical device
	Manufacturer
	Authorized representative
	CE mark reg. IVDD 98/79/EC
	Batch code
	Use by
	Cautions for safety and optimum product use
	Do not reuse

Endomedical™
1-800-282-6542
www.endomedical.ca

Manufactured for:
Endomedical Inc.
775 The Queensway Unit B
Toronto, Ontario
M8Z 1N1

Made in Korea

©2012 Endomedical. All Rights Reserved., Endomedical and Bravo are registered trademarks of Endomedical Inc.

 182100, 182101, 182102, 182104, 182105

 Manufacturer:
i-SENS, Inc. 465-6 Wolgye4-dong,
Nowon-gu, Seoul 139-845, Korea
www.i-sens.com

 **Medical Technology Promedt**
Consulting GmbH, Altenhofstrasse
80, D-66386 St. Ingbert, Germany

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©2012 Endomedical. Tous droits réservés., Endomedical et Bravo sont des marques de commerce déposées d’Endomedical Inc.